

ED | Inbound Patient – Set Arrival Time/Close Channel

01 Open Patient Channel

- From the *Patients* screen, tap the patient name or
- Scan the wristband

02 Set ED Arrival Time

- Scroll down to the *Emergency Department* card and tap on the *ED Arrival* time button.





03 Review or Assign Room (*Optional*)

- Tap *Edit* on the *Patient* card to add or update the *Room* number.

04 STOP Patient Channel

- When communication is no longer needed tap *STOP* to close the channel.


- 05 You can restart the channel, by scrolling to the bottom and tapping the *Restart Case* button.



ED | Inbound Patient – Set Arrival Time/Close Channel

01 Open Patient Channel

- From the *Patients* screen, tap the patient name or
- Scan the wristband.





ED | Inbound Patient – Set Arrival Time/Close Channel

01 Open Patient Channel

- From the *Patients* screen, tap the patient name or
- Scan the wristband.

02 Set ED Arrival Time

- Scroll down to the *Emergency Department* card and tap on the *ED Arrival* time button.



ED | Inbound Patient – Set Arrival Time/Close Channel



01 Open Patient Channel

- From the *Patients* screen, tap the patient name or
- Scan the wristband

02 Set ED Arrival Time

- Scroll down to the *Emergency Department* card and tap on the *ED Arrival* time button.

03 Review or Assign Room (*Optional*)

- Tap *Edit* on the *Patient* card to add or update the *Room* number.

David Miller
65 Years - Male

General
Transporting to
Grand View Regional

Patient DL Scan Edit

Statewide Wristband PMN FX A29517
Date of Service Jan 10, 2025 14:00
Name David Miller
Date of Birth Feb 08, 1960
ETA 14:30
Arrival Inbound EMS
Room -
Age 65 Years
Gender Male

Audio Images 1 ECGs 1

EMS

EMS Heartland EMS
Unit 09-5664521
First Contact 1

00:07:28 First Contact 00:00:00 ED Arrival Team Messages

David Miller - 65y M
Edit Patient

Patient DL Scan

Patient Movement Wristband#
PMN FX A29517

Name David Unknown
Miller Unknown

Date of Birth
Month Day Year
02 / 08 / 1960 Set Age

Age
65 Years Months Days

Gender
Female Male Unknown

Audio Images ECGs

Room ER4

ED | Inbound Patient – Set Arrival Time/Close Channel



01 Open Patient Channel

- From the Patients Screen, tap the patient name or
- Scan the Wristband

02 Set ED Arrival Time

- Scroll down to the Emergency Department card and tap on the ED Arrival Time button.

03 Review or Assign Room (*Optional*)

- Tap Edit on the Patient card to add or update the **Room** number.

04 STOP Patient Channel

- When communication is no longer needed, tap **STOP** to close the channel.

- 05 You can restart the channel by scrolling to the bottom and tapping the **Restart Case** button.

