EMS | Create Patient - Prehospital Notification



01) Scan Wristband - Create Patient

Tap **NEW** if not using Patient Movement wristband

02 Scan Driver's License - Demographics

- Enter demographics manually if not available
- Rapid: Unknown name, estimated Age, Gender

03 Photos

- Images: Patient, scene, paperwork...
- **ECG**: Photo or attach from monitor

⁰⁴ Patient Type

- Select General
- Upgrade as needed anytime

05 Voice to Text additional details

- Dictate your normal prehospital report
- "Eat the mic" Hold microphone close to mouth

Create Patient | Alert Hospital

- Of the Create: Tap Create to create channel now but alert hospital later
- 07 Alert Hospital Now: Tap Select Destination then Alert
 - Create button is replaced with Alert after destination is set

Note: All patient information can be added at any time, before or after patient creation.

03







